

# K.LASER

## GRANT YOURSELF THE BEST

 **BLUE**<sup>TM</sup>  
DENTAL 



# THE FIRST dental BLUE LASER both in clinics and aesthetics operating at 3 wavelengths

## THE FIRST LASER IN THE WORLD IN TERMS OF NUMBER OF APPLICATIONS IN THE DENTAL FIELD


**K-Laser Blue Dental** uses **445 nm** as its main wavelength which, unlike all other infrared lasers, is based on haemoglobin absorption and not on water absorption. This has many surgical advantages. Thanks to higher energy levels compared to infrared lasers, it has an increased germicidal power

The wavelengths, **970 nm** and **660 nm** support both superficial and deep tissue biostimulation.

**ADVANTAGES** No maintenance, 5-year warranty, low consumables cost and maximum safety. All the components can be sterilised.

### K-Laser Blue:

**445 nm - 660 nm - 970 nm**



### High surgical and germicidal activity

The 445nm wavelength interacts with the molecules at the systemic level. It is more easily absorbed by both haemoglobin and melanin, usually bigger than infrared laser. Blue laser is also effective on several bacteria strains, thus exploiting its antiseptic properties. It allows obtaining excellent results in terms of tissue incision both for surgical or vascular purposes.

### Faster tissue healing process

K-Laser Blue Dental wavelengths transfer the right amount of energy on the tissue surface. The energy is stored at a cellular level; this leads to an increase in both the proliferation and the metabolism. Thus, laser therapy allows obtaining excellent results in the tissue healing process.

### Cellular biostimulation and post-operative pain relief


Blood, the main transport system, provides the cells with all the elements that are necessary to their metabolism, such as oxygen and nutritional substances, thus removing catabolic products. The water contained in our body is mostly absorbed by the 970nm wavelength and a large part of the laser is transformed into heat. Thus, the deeper tissue layers become localised heat points, able to generate temperature gradients at a cellular level and to stimulate local microcirculation thanks to the increased cellular oxygen supply.

# K-LASER BLUE DENTAL

## New source of energy

**K-Laser Blue:** within the electromagnetic spectrum and among non-ionizing radiations, the most energetic wavelength is the one able to produce the maximum benefit for the tissues.

### THE ELECTROMAGNETIC SPECTRUM



# K-LASER BLUE DENTAL

## Excellency in surgery


---

**K-Laser Blue** has an interface that allows choosing among several treatments, ensuring the best surgical interventions and precise and unrivalled cutting performances.


**K-Laser Blue**, compared to other dental devices using only one infrared wavelength, uses blue light with a higher surgical performance, thus reducing thermal damage and offering optimal interaction with haemoglobin. Thanks to these features, it can guarantee a far better cutting precision compared to infrared laser.

Thanks to the non-absorption of water, surrounding tissues overheating is drastically reduced.

**K-Laser Blue** ensures the highest safety as the dental practitioner can avoid using scalpels and stitches thanks to the immediate haemostasis. The use of special sterilisable fibres avoids possible contamination and makes for a clean operating field.




**Infrared  
Laser**



# K-LASER BLUE DENTAL

## Applications

The way in which **K-Laser Blue** special wavelengths interact with the tissues as well as with the transferred energy determines the therapeutic and surgical indications below.



# K-LASER BLUE DENTAL

## Clinical tests


---

**K-Laser Blue Dental** can be used to treat infective (viral, bacterial and fungal) and traumatic wounds. In this case, the laser can be used to improve and reduce healing time after oral surgeries (removal of wisdom teeth and gingivectomies). Besides, it can also be used in symptomatic lesions with immune-mediated aetiology (recurrent mouth ulcers, lichen planus, pemphigus and pemphigoid). Lastly, the biostimulation and the analgesia deriving from the use of such wavelengths can be used in ATM dysfunctional pain syndromes.

### Lower lip herpes wound



Before




After two laser therapy sessions

### Mucositis of the buccal mucosa and the lingual margin



Before



After four laser therapy sessions

# K-LASER BLUE DENTAL

## Clinical tests

---

### Gingivectomy



Before the operation



End of operation



Check after two days

### Gingivectomy



Before the operation



End of operation



Check after 15 days

### Epulis



Before the operation



End of operation




Check after 4 days

# K-LASER BLUE DENTAL

## Clinical tests

---


### Fibroma



Before the operation



Treatment



End of operation



Check after 15 days


### Fibroma




Before the operation



End of operation



Check after 15 days




Check after 1 month

### Angioma




Before the operation



End of operation



Check after 15 days



Check after 1 month



# K-LASER BLUE DENTAL

## Aesthetic applications

---



Before



After



Before



After



Before




After




Before



After



Before



After



Before



After

# K-LASER BLUE DENTAL

K-Laser **intuitive software** is quick and precise and easily guides the dentist through each available treatment.

The parameters that determine the stage of each treatment, such as wavelength, energy, time, power and frequency, are automatically handled while leaving the operator the freedom to create new personalised protocols.




# K-LASER BLUE DENTAL

## Safe investment!

**K-Laser Blue Dental** is different from other laser equipment on the medical market for the wide choice of applications, the lack of maintenance and fast return on investments.

### Infrared laser devices



# K-LASER BLUE DENTAL

**K-Laser** has chosen to use semiconductors for their advantages compared to the other technologies on the market. Thus, **K-Laser**, proves to be the brand of the most reliable and compact devices in the world.

|   | K-Laser Blue              | Diode Laser   | Nd:YAG laser | Er:YAG laser      | Co2               |
|---|---------------------------|---------------|--------------|-------------------|-------------------|
| <b>Wavelength</b>                             | 445 nm + 660 nm<br>970 nm | 808 to 980 nm | 1064 nm      | 2940 nm           | 10600 nm          |
| <b>Model type</b>                             | Compact                   | Compact       | Large        | Extra large model | Extra large model |
| Maintenance costs                             | ▲                         | ▲             | △            | —                 | —                 |
| <b>Instructions</b>                           |                           |               |              |                   |                   |
| Surgery (soft tissues)                        | ▲                         | △             | △            | △                 | △                 |
| Vascular applications                         | ▲                         | —             | —            | —                 | —                 |
| Haemosthesis                                  | ▲                         | △             | —            | —                 | —                 |
| Periodontology                                | ▲                         | ▲             | △            | —                 | —                 |
| Endodontics                                   | ▲                         | ▲             | △            | ▲                 | —                 |
| Desensitisation                               | ▲                         | △             | △            | —                 | —                 |
| Biostimulation (LLLT)                         | ▲                         | ▲             | △            | —                 | —                 |
| Viral, bacterial, fungal infections           | ▲                         | △             | △            | —                 | —                 |
| Teeth whitening                               | ▲                         | △             | —            | —                 | —                 |
| Sealing                                       | ▲                         | —             | —            | —                 | —                 |
| Aesthetic applications                        | ▲                         | —             | —            | —                 | △                 |
| ▲ = Fit      △ = Partially fit      — = Unfit |                           |               |              |                   |                   |



# K-LASER BLUE DENTAL

## Technical file

**K-Laser Blue Dental** has been designed in compliance with directives 93/42/CEE and 2007/47/CE about medical products.

According to the norms in force, the device is classified as follows:

|                      |  |
|----------------------|--|
| Source type          | Integrated semiconductor laser module                                    |
| Laser system         | Class IV (compliant with IEC 60825-1)                                    |
| Product class        | Class IIb (compliant with directive 93/42/EEC)                           |
| Waivelength (nm)     | 445 nm $\pm$ 5nm; 660 nm $\pm$ 10 nm; 970 nm $\pm$ 15 nm                 |
| Overall power (W) CW | 7  |
| Emission Mode        | CW (continuous emission), pulsed, modality ISB, modality DISB            |
| Insulation Class     | Class II, type B (compliant with CEI EN 60601-1)                         |
| Beam Steering        | 660 nm $\pm$ 10 nm, max. 1mW   |
| DNRO                 | 12.46 m max  |
| Laser Activation     | wireless pedal   |
| Power supply         | Rechargeable battery and external power supply 100 - 240 VAC, 47 - 63 Hz |
| Display              | Full color, LCD touchscreen  |
| Size (W x L x H)     | 180 x 200 x 200 mm   |
| Body structure       | Aluminium body structure with special treatment                          |
| Handpiece            | Interchangeable in special metal. They can be sterilised                 |
| Fibers               | 200 $\mu$ m 320 $\mu$ m Multiuse sterilizable                            |
| Weight               | 2500 g (including the handpiece and the rechargeable battery)            |



# ACCESSORIES

## K-LASER BLUE DENTAL

---

Revolutionary innovation for advanced teeth laser surgery **K-Laser**



Optic fibre 200µm Art. **6255678**  
Optic fibre 320µm Art. **6255629**  
Sterilizable



Sterilizable tip  
Art. **MP419**



Handpiece body for contact surgery  
Art. **MP384**



Sterilizable sleeve for contact surgery  
Art. **MP385**

# ACCESSORIES

## K-LASER BLUE DENTAL

---



Handpiece body for non-contact surgery  
Art. MP386



Non-contact surgical spacers set  
Art. MP387A-MP387B-MP387C



Optic handpiece for therapy  
Art. MP383



ORL optics for therapy, sealing and whitening  
Art. MP388



Trolley with support for Tips  
Art. IM023



Wi-Fi Pedal  
Art. PF067



K-Laser suitcase  
Art. IM017A (optional)



Protective goggles for dentist and for patient  
Art. PF022P



Protective goggles with glass for patient  
Art. 6541523



Trolley with support for Tips  
Art. PF095 (optional)



For further information  
please visit the website  
[www.klaser.eu](http://www.klaser.eu)  
or scan this code to find  
out about the K-Laser  
Technology



K-Laser Blue Dental IT Rev. 1  
The company may make changes which it deems necessary without any notice. Reproduction, even only partial, is strictly forbidden

**KLASER**

ELTECH S.r.l.  
Strada Castagnole, 20/H  
31100 TREVISO - ITALY  
Tel. +39 0422 210430  
Fax +39 0422 297137  
com@klaser.it  
[www.klaser.eu](http://www.klaser.eu)

**KLASER**

**BLUE  
DENTAL II**

K-LASER and BLUE DENTAL  
are registered trademarks  
by Eltech srl



UNI ISO 9001:2008  
UNI ISO 13485:2012



0476